SOP:ATAC-seq for cryopreserved nuclei preparations from fresh tissueDate Modified:02/11/2020Modified by:Michelle Halstead (University of California, Davis)

The following protocol describes the isolation of nuclei from freshly harvested tissues, and subsequent preparation of ATAC-seq libraries for analysis of genome-wide open chromatin.

Chemicals Ordering Information

Item Dimethyl Sulfoxide (DMSO Hybri-Max 5 x 10 mL) D-Sucrose MgCl ₂ 1M (100 mL) Tris-HCl 1M pH 7.5 (1 L) Complete EDTA Free Protease inhibitor Tablets, mini Phosphate buffered saline Tris-HCl NaCl IGEPAL CA-630 DPBS, no calcium, no magnesium Nextera DNA Sample Prep Kit SsoFast™ EvaGreen® Supermix Qiagen MinElute PCR Purification Kit PCR Primers (custom sequences: see Table 1 below)	Catalog Number D2650 BP220-1 AM9530G 46-030-CM 04-693-132-001 10010023 10812846001 S9888 I8896 14190250 FC-121-1030 1725200 28004	Manufacturer Sigma-Aldrich Fisher Scientific Ambion Mediatech, Inc. Roche Applied Science ThermoFisher Sigma-Aldrich Sigma-Aldrich Sigma-Aldrich ThermoFisher Illumina Bio-Rad Qiagen IDT
Qiagen MinElute PCR Purification Kit PCR Primers (custom sequences; see Table 1 below) Milli-Q or Molecular Biology Grade Sterile Water	28004	Qiagen IDT

Table 1. Oligo designs. A list of ATAC-seq oligos used for PCR. Primers were ordered from IDT. Stock solutions were prepared according to instructions from IDT (100 μ M) and were stored at -20°C. Working solutions (25 μ M) were also stored at -20°C.

Primer	Sequence (5' to 3')
1	AATGATACGGCGACCACCGAGATCTACACTCGTCGGCAGCGTCAGATGTG
2A	CAAGCAGAAGACGGCATACGAGATCTAGTACGGTCTCGTGGGCTCCGGAGATGT
2B	CAAGCAGAAGACGGCATACGAGATTTCTGCCTGTCTCGTGGGCTCGGAGATGT
2C	CAAGCAGAAGACGGCATACGAGATGCTCAGGAGTCTCGTGGGCTCGGAGATGT
2D	CAAGCAGAAGACGGCATACGAGATAGGAGTCCGTCTCGTGGGCTCGGAGATGT
2E	CAAGCAGAAGACGGCATACGAGATCATGCCTAGTCTCGTGGGCTCGGAGATGT
2F	CAAGCAGAAGACGGCATACGAGATGTAGAGAGGTCTCGTGGGCTCGGAGATGT
2G	CAAGCAGAAGACGGCATACGAGATCCTCTCTGGTCTCGTGGGCTCGGAGATGT
2H	CAAGCAGAAGACGGCATACGAGATAGCGTAGCGTCTCGTGGGCTCGGAGATGT
21	CAAGCAGAAGACGGCATACGAGATCAGCCTCGGTCTCGTGGGCTCGGAGATGT
2J	CAAGCAGAAGACGGCATACGAGATTGCCTCTTGTCTCGTGGGCTCGGAGATGT
2К	CAAGCAGAAGACGGCATACGAGATTCCTCTACGTCTCGTGGGCTCGGAGATGT
2L	CAAGCAGAAGACGGCATACGAGATATCACGACGTCTCGTGGGCTCGGAGATGT
2M	CAAGCAGAAGACGGCATACGAGATACAGTGGTGTCTCCGTGGGCTCGGAGATGT
2N	CAAGCAGAAGACGGCATACGAGATCAGATCCAGTCTCGTGGGCTCGGAGATGT
20	CAAGCAGAAGACGGCATACGAGATACAAACGGGTCTCGTGGGCTCGGAGATGT
2P	CAAGCAGAAGACGGCATACGAGATACCCAGCAGTCTCGTGGGCTCGGAGATGT
2Q	CAAGCAGAAGACGGCATACGAGATAACCCCTCGTCGTGGGGCTCGGAGATGT
2R	CAAGCAGAAGACGGCATACGAGATCCCCAACCTGTCTCGTGGGCTCGGAGATGT
2S	CAAGCAGAAGACGGCATACGAGATCACCACACGTCTCGTGGGCTCGGAGATGT
2T	CAAGCAGAAGACGGCATACGAGATGAAACCCAGTCTCGTGGGCTCGGAGATGT
2U	CAAGCAGAAGACGGCATACGAGATTGTGACCAGTCTCGTGGGCTCGGAGATGT
2V	CAAGCAGAAGACGGCATACGAGATAGGGTCAAGTCTCGTGGGCTCGGAGATGT
2W	CAAGCAGAAGACGGCATACGAGATAGGAGTGGGTCTCGTGGGCTCGGAGATGT
2X	CAAGCAGAAGACGGCATACGAGATTCGCCTTAGTCTCGTGGGCTCGGAGATGT

Materials List

Dounce 7mL Tissue Grinder with PYREX Pestles, Corning (VWR Cat# 22877-280) GentleMACS dissociator from Mitenyi Biotec. GentleMACS C-tubes (Mitenyi Biotec Cat# 130-093-237) 100 µm Steriflip 50mL Disposable Vacuum Filter System (Millipore Cat# SCNY00100) Nalgene Cryo 1°C Freezing Container (Cat# 5100-0001) Liquid Nitrogen Storage Graduated pipets (5, 10, 25, 50 mL) Micropipettes and tips (P2, P20, P200, P1000) 500 mL Corning 0.2 µM Filter System (Cat# 430758) 1 L Corning 0.2 µM Filter System (Cat# 430186) 0.2-ml PCR tubes 1.5-ml Eppendorf tubes 15 mL Corning Polypropylene Conical Centrifuge Tubes (Cat# 430766) 50 mL Corning Polypropylene Conical Centrifuge Tubes (Cat# 430828) 20 µM Steriflip 50mL Disposable Vacuum Filter System (Millipore Cat# SCNY00020) CryoTube Vials, 1.8 mL (Nunc Cat# 368632) Microscope (preferably with phase contrast) Hemacytometer (Hausser Scientific Cat# 1492) Refrigerated centrifuge (swing-out rotor) compatible with 15 mL Conical and 1.5-ml Eppendorf tubes Thermal Mixer (ThermoFisher Cat# 13687711) PCR thermal cycler 2100 Bioanalyzer Instrument from Agilent Bioanalyzer High Sensitivity DNA Kit (Agilent Cat# 5067-4626) PippinHT size-selection system from Sage Science PippinHT 3% Agarose DNA size-selection cassette, 100-250bp (Sage Science Cat# HTG3010)

Stock Reagents

Unless otherwise noted, all buffers and stock solutions should be pre-chilled to 4°C (on ice) prior to use.

Sucrose buffer				
Final concentration	Stock concentration	Amount used from stock		
250mM D-Sucrose	0.5M D-Sucrose	250 mL		
10mM Tris-HCl, pH 7.5	1M Tris-HCl, pH 7.5	5 mL		
1mM MgCl ₂	1M MgCl ₂	0.5 mL		
Molecular Biology Grade sterile H ₂ O to 500 mL				
Filter sterilize with 500 mL 0.2 μ M Filter System. Store at 4°C. Add Complete Protease Inhibitor				
Tablets (1 per 50mL solution) just prior to use.				

ATAC-seq cell lysis buffer

Final concentration	Stock concentration	Amount used from stock	
10mM Tris-HCl, pH 7.4	1M Tris-HCl, pH 7.4	50 μL	
10mM NaCl	1M NaCl	50 μL	
3mM MgCl ₂	1M MgCl ₂	15 μL	
0.1% (v/v) IGEPAL CA-630		5 μL	
Molecular Biology Grade sterile H ₂ O to 5 mL			
Good for up to 1 week; recommend making fresh each time if possible.			

Prior to nuclei isolation:

- 1. Add protease inhibitor tablets to sucrose Buffer (1 tablet per 50 mL solution) and solubilize. Keep on ice.
- 2. Pre-label and open all sterile filters and GentleMACS C-tubes.
- 3. Work quickly using reagents maintained at appropriate temperatures.

Cryopreservation of nuclear prep samples using gentleMACS dissociator

- Collect ~3mm square piece or roughly 1 gram of fresh tissue <u>Note</u>: Less material could be acceptable for ATAC-seq, which should work with a minimum of about 500 cells.
- 2. Transfer minced tissue into a gentleMACS C-tube with 10 ml of ice-cold sucrose buffer.
- 3. Mince tissue with scissors in in the gentleMACS C-tube. Pieces should be cut up to about 1/10 of the size of the square.
- 4. Homogenize tissue using gentleMACS Dissociator Program "E.01c Tube" *twice*.
- 5. Filter homogenate using 100 μm Steriflip Vacuum Filter System. <u>Note:</u> Pour slowly and check for large pieces as they can clog the filter. Gently scraping the filter surface can help dislodge clogs.
- 6. Bring volume to 9.9 ml with Sucrose buffer.
- 7. Add 1.1 mL DMSO to samples (10% final concentration), pipetting several times to adequately mix.
- 8. Aliquot into 5 cryotube vials, with approximately 1.5 mL in each vial.
- 9. Freeze at -80°C overnight in Nalgene Cryo 1°C Freezing Container, then move to -80°C freezer or -135°C liquid nitrogen for long-term storage.

Cryopreservation of nuclear prep samples using Dounce homogenizer

- 1. Ideally, keep equal weight between tissues (~3mm square tissue piece or roughly 1 gram)
- 2. Mince tissue with razor blade or scissors in polystyrene weighing dish.
- 3. Homogenize tissue with 3 mL Sucrose Buffer per gram tissue in Dounce tissue grinder.
- 4. Dounce approximately 5 times, slowly and smoothly, with loose (A) pestle.
- 5. Filter homogenate using 100 μm Steriflip Vacuum Filter system.
- 6. Bring up to 2.7 mL with Sucrose Buffer.
- 7. Add 0.3 mL DMSO to samples (10% final concentration), pipetting several times to adequately mix.
- 8. Aliquot into cryotube vials, freeze at -80°C overnight in Nalgene Cryo 1°C Freezing Container, then move to -80°C freezer or -135°C liquid nitrogen for long-term storage.

Prior to ATAC-seq:

- 1. Chill centrifuge *with swing-out rotor* to 4°C.
- 2. Prepare and chill ATAC-seq cell lysis buffer to 4°C. Ideally, cell lysis buffer should be prepared fresh.
- 3. Thaw cryopreserved nuclei preparations on ice.

ATAC-seq library generation from cryopreserved nuclei preparations

- 1. After thawing on ice, transfer nuclei preparations to 15 mL conical centrifuge tubes and centrifuge for 5 min, 500 RCF, 4°C.
- 2. Aspirate supernatant, gently resuspend pellet in 1 mL ice-cold DPBS, and centrifuge for 5 min, 500 RCF, 4°C.
- 3. Aspirate supernatant, gently resuspend pellet in 1 mL ice-cold ATAC-seq cell lysis buffer, and centrifuge for 10 min, 500 RCF, 4°C.
- Aspirate supernatant and gently resuspend pellet in 1 mL cold DPBS. <u>Note:</u> Depending on the abundance of cells in the preparation, resuspending in smaller or greater volumes of DPBS may be necessary to accurately count nuclei.
- 5. Determine nuclei concentration using a hemocytometer.

<u>Note:</u> For some samples, cells will be more commonly observed than nuclei. This may indicate incomplete cell lysis, which can be checked by staining with propidium iodide and Hoechst's.

- 6. Aliquot between 50,000 and 500,000 cells to 1.5 mL Eppendorf tubes and centrifuge for 5 min, 500 RCF, 4°C.
- 7. Carefully aspirate supernatant and gently resuspend pellet in ice-cold 50 μ L transposition mix (25 μ L TD buffer from Nextera DNA Sample Prep Kit, 2.5 μ L TDE1 enzyme from Nextera kit, 22.5 μ L ddH₂O).
- 8. Incubate nuclear pellet with transposition mix for 60 min, 37°C, 300 rpm.
- 9. Purify transposed DNA with MinElute PCR purification kit, and elute DNA in 10 µL Buffer EB.

Optional stopping point: Store transposed DNA at -20°C.

10. Add 40 μL PCR master mix (25.4 μL SsoFast[™] EvaGreen[®] Supermix, 13 μL ddH₂O, 0.8 μL 25 μM Primer 1, 0.8 μL 25 μM Primer 2) to 10 μL eluted DNA and cycle as follows:

1x	5 min 30 sec	72°C 98 °C
10 - 13x	10 sec 30 sec 1 min	98 °C 63 °C 72 °C

<u>Note</u>: Fewer PCR cycles is desirable to minimize PCR duplication. Furthermore, if sufficient DNA concentration is obtained from fewer PCR cycles, it is a good indicator of library complexity. However, more cycles may be necessary to obtain sufficient DNA concentration in some samples.

- 11. Purify libraries using MinElute PCR purification kit, and elute DNA in 10 μL Buffer EB.
- 12. Quantify and characterize libraries using Agilent Bioanalyzer High Sensitivity DNA chip.

<u>Note:</u> Libraries with concentrations less than 6 nM (150-250 bp range) should be considered for additional cycles of PCR. Keep in mind, this will increase PCR duplication rate after sequencing. If possible, repeating library construction may be preferable.



- Size-select libraries for subnucleosomal fragments (150-250 bp) on the PippinHT system using a 3% Agarose DNA size-selection cassette.
- 14. Quantify library concentration and verify size-selection on Agilent Bioanalyzer High Sensitivity DNA chip.



15. Pool libraries with compatible primer indexes and submit for sequencing using a platform such as NextSeq 500 to generate paired-end reads. Ideally, 40 million reads should be obtained per each library in the pool (not including low-quality, mitochondrial, or duplicate alignments).