

RNA extraction from Atlantic Salmon head kidney tissue

Product	Company	Reference
Water Molecular Biology Grade (1000ml)	Thermo Scientific	SH30538.03
Tungsten Carbide Beads	Qiagen	69997
RNaseZap®, 250 mL, s, cat. no. , £47.20	Life technologies	AM9780
Sterile Microtubes	Axygen	MCT-150-C-S
TRI Reagent®	Sigma	T9424-200ML
RNAlater®	sigma	R0901-500ML
Chloroform	sigma	C2432-1L
2-Propanol	sigma	I9516-500ML
Molecular grade 100% Ethanol	sigma	2483

- Thaw head kidney in RNAlater kept at -80°C
- Take a small piece of the tissue (10-20mg), get rid of excess RNAlater using a clean paper towel.
- Add the tissue into 1.5 mL of cold TRIzol in a 2 mL tubes with 2 Tungsten (3 mm) beads
- Homogenise at frequency 30.0 l/sec for 2.5 min, 30 sec breaks
- Incubate on ice for 5 minutes to permit complete dissociation of the nucleoproteins complex
- Add 0.3ml of chloroform, vortex to mix well
- Incubate for 5 minutes on ice
- Centrifuge the sample for 15 minutes at 12,000xg at 4°C
 - The mixture separates into a lower red phenol-chloroform, and interphase, and a colorless upper aqueous phase
- Transfer the aqueous phase containing the RNA to a new tube, be careful to not take any of the interphase layer
- Add 0.7ml of isopropanol to the aqueous phase, mix well by vortexing
- Incubate for 10 minutes and centrifuge for 10 minutes at 12000xg at 4°C
 - Total RNA precipitate forms a white gel-like pellet at the bottom of the tube
- Discard the supernatant
- Resuspend the pellet in 1ml of 80% ethanol per 1ml of trizol used for lysis, vortex until pellet detached
- Centrifuge for 5 minutes at 7500xg at 4°C
- Discard the supernatant and repeat the washing with 80% ethanol two more times

- Pour away alcohol (watch the pellet), then use 20 μ L pipette to remove alcohol (do not touch the pellet), close tubes & spin down till 3500xg. Use 10 μ L tips to get rest of alcohol (do not touch the pellet), close the tubes, spin down till 3500xg, then use 10 μ L tips again for the final alcohol removal (make sure that all alcohol is gone, otherwise the alcohol in the sample may inhibit qPCR), close tubes.
- Let the RNA pellet air dry for 10 minutes on ice, or if needed, put tubes on hot plate 65°C for 1 min to make sure all alcohol evaporates
- Dissolve the pellet in 100ul of nuclease free water
- Flip the tube with the finger and leave the tubes for 30 min on wet ice
- If pellet is not dissolved after 30 minutes of incubation, incubate the tubes at 65C on hot plate for 10 minutes to properly dissolve the pellet
- Transfer the RNA solution into a new RNAase free tube and check quantity and quality on nanodrop and bioanalyzer
- Store on -80C for longer storage