

UNIVERSITÀ
DEGLI STUDI
DI PADOVA

ATAC-seq from seabass fresh cells

Overview: This protocol describes the method used to prepare ATAC-seq libraries from seabass cell cultures prepared as described in

https://data.faaang.org/api/fire_api/samples/UNIPD_SOP_Leucocytes_Isolation_and_Stimulation%20Dlab_20200518.pdf

Day 1 – Cells collection, transposition and clean up

Cells collection

NB: In order to reflect seabass cell osmolarity, PBS 1X was adjusted to 350 mOsm/kg by adding 2M NaCl.

1. Harvest all cells by softly pipetting the cell suspension with a WO tip and then collect them in a 2 ml Eppendorf tube.
2. Pellet the cells by centrifugation at 610g for 5 min at 4°C.
3. Eliminate the supernatant and resuspend the pellet with 1mL of cold **PBS 1X (350 mOsm/kg)**
4. Centrifuge the suspension at 610g for 5 min at 4°C
5. Eliminate the supernatant and resuspend the pellet with 500 ul of cold **PBS 1X (350 mOsm/kg)**
6. Count the cells with hemocytometer (5ul diluted cell suspension + 5 ul PBS + 10 ul Trypan blue)

Nuclei isolation, transposase reaction and clean-up

7. Put a volume of cell suspension corresponding to 50.000 cells (from step 6) in a 1.5 ml vial
8. Pellet 50,000 cells by centrifugation at 610g for 5 min at 4 °C
9. Add 50 ul cold **ATAC-RSB-L** buffer and incubate on ice for 3 minutes
10. Wash out lysis with 1 mL cold **ATAC-RSB-W**, invert the tube at least 3 times to mix
11. Centrifuge nuclei at 800 g for 10 minutes at 4 °C
12. Aspirate supernatant carefully using a p1000 followed by p100 pipette. The nuclei pellet is usually very small and can only be observed as an opaque “shadow”. Take care not to also aspirate the nuclei pellet

13. Make the following transposase reaction mix:

Reagent	Volume per sample (ul)
2X Tagment DNA (TD) Buffer	25
Transposase	1.25
PBS 1X	16.5
Digitonin (1 %)	0.5
Tween-20 (10 %)	0.5
Nuclease free H ₂ O	6.25
Total	50

14. Resuspend the pellet in 50 ul transposase reaction mix by pipetting up and down
15. Incubate the reaction at 37 °C for 30 minutes in a thermomixer with 1000 RPM mixing
16. Stop the reaction and purify the resulting DNA fragments with MinElute PCR purification kit
16.1 Elute DNA in 23 ul EB buffer
17. **SAFE STOP POINT** - Store the eluted DNA at -20°C.

Day 2 – Library amplification

Initial PCR amplification

18. Save the following program (ATAC-PRE) on a thermal cycler with a heated lid:
72°C for 5'
98°C for 30''
98°C for 10''
63°C for 30''
72°C for 1'
4°C forever
} 5 cycles
19. Set up the following PCR reaction:
19.1. Add 25ul of NEBNext® Ultra™ II Q5® Master Mix to each tube.
19.2. Add 5 ul of Index from IDT® for Illumina Nextera DNA Unique Dual Indexes (UDI) Set
19.3. Add 20 ul of transposed DNA (from step 16)

Reagent	Vol per sample (ul)
IDT® for Illumina Nextera DNA UDI	5
NEBNext Ultra II Q5 Master Mix	25
Transposed DNA	20
Total	50

20. Mix reagents, close tubes and centrifuge at 280 g for 1 min
21. Place the plate on the preprogrammed thermal cycler and run ATAC-PRE

qPCR assessment of final PCR cycles

22. Using 5 ul (10 %) of the Initially-amplified product, assemble 15 ul qPCR reactions in a qPCR plate to determine the appropriate number of additional cycles needed.

22.1. Prepare a qPCR mix of nuclease free water, NEBNext® Ultra™ II Q5® Master Mix and 25X SYBRGreen for n samples +2 (n samples + 1 NTC + 1 additional).

Reagent	Vol per sample (ul)
Nuclease free water	3.76
NEBNext Ultra II Q5 Master Mix	5
25X SYBRGreen (diluted*)	0.4
Total	9

* SYBR green is provided at 10,000 X. Make 1ul aliquots and freeze, then on the day of use add 399ul H2O and mix well.

- 22.2 Distribute 9 ul of qPCR mix in each well.
 22.3 Add 1 ul of Nextera UDI index to each sample. IMPORTANT, this must be the same as was used for Initial PCR reaction.
 22.4 Add 5 ul of Initially amplified DNA (from step 21) or 5 ul of water (for NTC)
23. Seal the plate and centrifuge at 280g for 1 min
 24. Run the following program on a Real-time thermal cycler:

98°C for 30''
 98°C for 10''
 63°C for 30'' } 20 cycles
 72°C for 1'
 30°C – Stop

25. Determine the required number of additional cycles to amplify. The number of cycles should equal $\frac{1}{4}$ of max fluorescence (Figure 1). This is to avoid PCR bias.




Figure 1. qPCR amplification plot showing $\frac{1}{4}$ of max fluorescence = 3 PCR cycles

26. Save the following program (ATAC-FINAL) on a thermal cycler with a heated lid

98°C for 30''
 98°C for 10''
 63°C for 30'' } n° cycles determined in step 25
 72°C for 1'
 4°C forever

27. Short-spin the PCR tubes containing your Initially-amplified product (45ul).

28. Place the plate on the thermal cycler and run ATAC-FINAL.

SAFE STOP POINT - After Final PCR amplification, keep PCR product at 4°C until you proceed with size selection.

Same day or another day – Library size selection

Library size selection

29. Resuspend Ampure XP beads by vortexing well (> 1 min)
30. Add 0.55x volume (24.75 ul) of beads to sample. Mix well by pipetting. Be gentle to avoid bubbles
31. Incubate at room temperature for 10 min
32. Place on magnetic rack and let it stand for 5 min (until the supernatant is clear of beads)
33. **Transfer the supernatant** to a new tube/plate and add another 1.3x original volume (58.5 ul) of Ampure beads to the supernatant
34. Mix well by pipetting. Be gentle to avoid bubbles
35. Incubate at room temperature for 10 min
36. Place on magnetic rack and let it stand for 5 min (until the supernatant is clear of beads)
37. Remove and **discard supernatant**
38. Wash the beads 2x with freshly made 80 % ethanol while plate/tube is still on the magnetic rack:
 - 46.1 Add 100 ul 80 % ethanol over beads
 - 46.2 Wait 1 minute
 - 46.3 Remove ethanol
 - 46.4 Repeat 1x
39. Remove samples from the magnet and allow the tubes to air dry (30sec-2min)
40. Add 16.5 room temperature TET Buffer. Resuspend beads by pipetting
41. Rehydrate at room temperature for a minimum of 2 minutes
42. Place on magnetic rack and let it stand for 5 min (until the supernatant is clear of beads)
43. Transfer the supernatant (eluted DNA) to a LoBind Eppendorf tube or 96 well plate

Library quality assessment

44. Measure DNA library concentration with Qubit High Sensitivity Kit
45. Use 1 ul of diluted library to validate DNA fragment size distribution with Bioanalyzer High Sensitivity DNA kit .

The DNA fragment size distribution should to some extent follow a nucleosome pattern (Figure 2) with the most prominent peak being at about 200 bp.




Figure 2. Example of Bioanalyzer trace after library preparation and size selection

BUFFERS

Buffers that can be prepared in advance (store at 4°C)

ATAC-RSB (store at 4°C)	Final conc.	Vol for 50 mL
1 M Tris pH 7.4	10 mM	500 µl
5 M NaCl	10 mM	100 µl
1 M MgCl ₂	3 mM	150 µl
dH ₂ O		49.25 mL

Detergents		
Digitonin is supplied at 2% in DMSO. Dilute 1:1 with water to make a 1% (100x) stock solution. Avoid more than 5 freeze thaw cycles. Can be kept at -20°C for up to 6 months.		
Tween-20 is supplied at 10%. Use at this concentration (100x stock). Store at 4°C		
NP40 is supplied at 10%. Use at this concentration (100x stock). Store at 4°C		

TET buffer (16.5 µl per sample) (store at 4°C)	Final conc.	Vol per sample (µl)	Vol for 100 samples
Tris-HCL pH 8.0 (1 M)	10 mM	0.165	16.5
EDTA (0.5 M)	1 mM	0.033	3.3
Tween-20 (10 %)	0.05%	0.0825	8.25
dH ₂ O		16.2195	1621.95

Buffers to be prepared the same day (keep on ice)

ATAC_RSB_L	Final conc.	Vol per sample (µl)
ATAC-RSB	1x	48.50
1% Digitonin	0.01%	0.50
10% Tween-20	0.10%	0.50
10% NP-40	0.10%	0.50
TOTAL		50.00

ATAC-RSB-W	Final conc.	Vol per sample (µl)
ATAC-RSB	1x	990.00
10% Tween-20	0.10%	10.00
TOTAL		1000.00