



Total RNA extraction from Atlantic salmon Head Kidney Leukocytes

Product	Company	Reference
Water Molecular Biology Grade (1000ml)	Thermo Scientific	SH30538.03
Tungsten Carbide Beads	Qiagen	69997
RNaseZap®, 250 mL, s, cat. no. , £47.20	Life technologies	AM9780
Sterile Microtubes	Axygen	MCT-150-C-S
TRI Reagent®	Sigma	T9424-200ML
RNAlater®	sigma	R0901-500ML
Chloroform	sigma	C2432-1L
2-Propanol	sigma	I9516-500ML
Molecular grade 100% Ethanol	sigma	2483

RNA extraction from 6/12 cell plates

- Centrifuge the plate at 500g for 5 minutes to make sure that cells are stuck on the bottom
- Put sterile papers on the sink and quickly invert the plate emptying the growth media without losing cells. Make sure you have removed all the media.
- In cell culture hood, add trizol depending on the cell count. Add 0.5ml of trizol into each well if you 1-2million cell count
- Pipette up and down all over the well to get all the cells attached to the bottom, put trizol and cell mixture in to a 1.5/2ml Eppendorf tube

You can freeze the samples at this point for later extraction, otherwise proceed with the protocol

- Incubate for 5 minutes to permit complete dissociation of the nucleoproteins complex
- Add 0.2ml of chloroform per 1ml of trizol used for lysis, vortex to mix well
- Incubate for 2-3 minutes
- Centrifuge the sample for 15 minutes at 12,000xg at 4C
 - The mixture separates into a lower red phenol-chloroform, and interphase, and a colorless upper aqueous phase
- Transfer the aqueous phase containing the RNA to a new tube, be careful to not take any of the interphase layer
- Add 0.5ml of isopropanol to the aqueous phase, per 1ml of trizol used for lysis, mix well by vortexing
- Incubate for 10 minutes and centrifuge for 10 minutes at 12000xg at 4C
 - Total RNA precipitate forms a white gel-like pellet at the bottom of the tube

- Discard the supernatant
- Resuspend the pellet in 1ml of 80% ethanol per 1ml of trizol used for lysis, vortex until pellet detached
- Centrifuge for 5 minutes at 7500xg at 4C
- Discard the supernatant and repeat the washing with 80% ethanol two more times
- Pour away alcohol (watch the pellet), then use 20 μ L pipette to remove alcohol (do not touch the pellet), close tubes & spin down till 3500xg. Use 10 μ L tips to get rest of alcohol (do not touch the pellet), close the tubes, spin down till 3500xg, then use 10 μ L tips again for the final alcohol removal (make sure that all alcohol is gone, otherwise the alcohol in the sample may inhibit qPCR), close tubes.
- Let the RNA pellet air dry for 10 minutes on ice, or if needed, put tubes on hot plate 65°C for 1 min to make sure all alcohol evaporates
- Dissolve the pellet in 40ul of nuclease free water
- Flip the tube with the finger and leave the tubes for 30 min on wet ice
- If pellet is not dissolved after 30 minutes of incubation, incubate the tubes at 65C on hot plate for 10 minutes to properly dissolve the pellet
- Transfer the RNA solution into a new RNAase free tube and check quantity and quality on nanodrop
- Store on -80C for longer storage