

Standard operating procedure to isolate genomic DNA from blastocysts

Research Institute for Farm Animal Biology (FBN)

Institute of Genome Biology, Genome Physiology Unit



RESEARCH INSTITUTE FOR
FARM ANIMAL BIOLOGY

DNA isolation

Samples

- frozen blastocysts

Equipment

- Centrifuge and benchtop centrifuge
- Qubit
- Pipets
- Stand
- Styrofoam box with ice
- Heating block

Plastics

- RNase/DNase-free tubes
- Pipet tips

Kits & Chemicals

- Quick-DNA Microprep Plus Kit chemicals
- Proteinase K from Quick-DNA kit.
- Qubit dsDNA HS Assay Kit

Preparation

- Pre-heat the heating block to 55 °C
- Dissolve the neophylised 5 mg Proteinase K in 260 µl Proteinase K Buffer from the kit to a concentration of 20 mg/ml
- Transfer some DNA Elution Buffer into a 1.5 ml tube for better preheating in the heating block

Procedure

1. Resuspend the blastocysts with 50 µl DNA Elution Buffer and transfer to a new 1.5 ml tube.
2. Add 50 µl BioFluid & Cell Buffer (Red) to the suspension. Also flush the 0.2 ml cell storage tube with the buffer.
3. Add 5 µl Proteinase K and vortex the tube for 10 to 15 seconds.
4. Incubate the mix at 55 °C in the heating block for 10 minutes.
5. Preheat the DNA Elution Buffer to 65 °C in the heating block.
6. Add one volume 105 µl Genomic Binding Buffer and vortex 10 to 15 seconds.
7. Centrifuge briefly.
8. Transfer the lysate to the Zymo Spin IC-XM Column and centrifuge at 12.000 g for 5 minutes at room temperature.
9. Discard the flowthrough and transfer the column in a new collection tube.

10. Add 200 µl DNA Pre-Wash Buffer to the column and centrifuge at 12.000 x g for 1 minute at room temperature.
11. Discard the flowthrough and transfer the column in a new collection tube.
12. Add 700 µl g-DNA Wash Buffer and centrifuge at 12.000 x g for 1 minute at room temperature.
13. Discard the flowthrough and transfer the column in a new collection tube.
14. Add 200 µl g-DNA Wash Buffer and centrifuge at 12.000 x g for 1 minute at room temperature.
15. Discard the flowthrough and transfer the column in a clean 1.5 ml tube.
16. Pipette 11 µl of 65 °C warm elution buffer directly onto the column matrix without touching the matrix.
17. Incubate for 5 minutes at room temperature.
18. Centrifuge at full speed for 1 minute.
19. Transfer the eluate one more time to the column matrix.
20. Incubate for 3 minutes at room temperature.
21. Centrifuge at full speed for 1 minute at room temperature.

For storage keep the samples at – 20 °C. Use the Qubit dsDNA HS Assay Kit to determine the DNA concentration with Qubit.