



Dechorionation of sea bass embryos and cleaning of chorion constituents

UNIPD 12-10-2020

<u>Overview:</u> This protocol describes a method used to extract sea bass embryos from eggs, and clean them from other egg constituents, such as proteins and oils which can interfere with downstream library preparation for ATAC-seq and ChIP-seq.

Consumables:

- Sterile 50 ml plastic tubes
- Watchmaker forceps
- 1.5ml Eppendorf tubes
- Pronase 1mg/ml in seabass-osmolar PBS 1X, (stock: 30 mg/ml)
- Seabass-osmolar PBS 1X (add 1100 μl of NaCl 2 M to 50 ml of PBS)
- Seabass-osmolar DMEM
- Seabass-osmolar FBS (heat inactivated)
- proteinase inhibitor cocktail tablet (Complete PIC tablet, EDTA free; Roche)
- 15ml Falcon tubes
- Rotator system (HulaMixer[™] Sample Mixer)
- Dissection microscope
- 70 µm filters
- 6-well cell culture plates (and/or Petri dish)
- Single concave microscope slides
- Cooled incubator
- Vortex
- Fixed angle centrifuge
- Crushed ice

Seabass osmolar Resuspension buffer (RSB)

5 ml FBS 5 ml DMEM 1 proteinase inhibitor cocktail tablet (Complete PIC tablet, EDTA free; Roche) 105 μl NaCl 2M

Seabass osmolar Cryopreservant solution

- 4 ml FBS 5 ml DMEM 1 DMSO 91 μl NaCl 2M
- 2.1 Place required number of eggs in a sterile 50ml plastic tubes with 6ml of pronase 1mg/ml in seabass-osmolar PBS 1X.
- 2.2 Rotate for 30 min at 20RPM at 20°C in a cooled incubator.
- 2.3 Eliminate pronase solution and add 6 ml of cold seabass-osmolar PBS 1X.
- 2.4 Transfer one, two or three eggs to a single concave microscope slide under the microscope with two or three cold PBS 1X drops.
- 2.5 Dechorionate using watchmaker forceps. Take care not to damage blastula/gastrula.
- 2.6 Use light mechanical pressure from the forceps to separate blastula from yolk.
- 2.7 Using a p200 set to 150ul, with a wide-bore tip, transfer blastula to a 1.5ml Eppendorf tube on ice.
- 2.8 Repeat for required number of eggs for each replicate, pipetting sample up and down.
- 2.9 Centrifuge on fixed angle centrifuge, 300rfc for 5min @ 4°C
- 2.10 Remove supernatant.

The embryos are ready for homogenization and cell isolation for ATACseq protocol from fresh tissue or for formalin fixing for ChIPseq protocol.